

FUNWEIGH™ (FUNNY TALKING SCALE) MANUAL



Congratulations on purchasing a genuine My Weigh Scale.

You now own a Personal Scale of proven technical design and highest quality.

We hope that you will enjoy the use of your new Personal Digital Scale.

You can visit us at www.Myweigh.com

Function

Measures body weight and with voice reading and funny commands.

Technical specifications:

Weighing capacity:	200kg / 440lb
Minimum weighing limit:	2kg / 4.4 lb
Weighing graduation:	0.1kg / 0.2 lb
Weighing deviation:	± (1%+1d)
Optimum temperature range:	20°C ± 15°C

Operation instructions:

1. Inserting batteries:
Insert 4 AA batteries into the battery compartment on the bottom (underside) of the scale. Be sure the polarity (+ / -) is correct when installing batteries
2. Taking measurement
 - a. Step onto the measurement platform. The scale will say "Hello". Now step off the platform and wait a few moments
 - b. Once the scale is stable, it'll say "It's ready" and show "0.0 lb" (or 0.0 kg). Then you may step onto the measurement platform again. Your weight will be displayed and spoken once the reading is steady, followed by a funny voice command.
 - c. You may step off of the platform. The scale will say "Good-bye" and turn off automatically after about 10 secs.
3. Power off
The power will shut down automatically after no operation is detected for about 10 seconds.
4. Switching on/off the talking function
Use the switch on the bottom of the platform to switch on/off the speaker function.
5. Unit conversion
Turn the Scale ON, "0.0" will show on the display" press the kg/lb button on the front of the scale (next to the LCD display) to select the desired unit (LB or KG)
6. Important notes
 - ***If the Scale ever malfunctions please try installing new batteries, sometimes a low battery or defective battery may cause inaccurate readings or other errors
 - When the batteries run very low, the warning icon will appear on the display. When this happens, be sure to replace the battery immediately on the underside of the scale.
 - When over-loaded, "FFFF" will be displayed.
 - Use the unit on hard, flat, dry and secure surface only.
 - Place the unit in an area free from high temperature or high humidity. To clean, use dampened cloth and mild detergent.
 - Avoid excessive impact or vibration to the unit.
 - If you do not intend to use this unit for a long period of time, remove the battery before storage.
 - Do not attempt to disassemble the unit, otherwise the warranty will become void.
 - The unit is intended for home and personal use only (not legal for commerce).
 - **This scale must only be used on a hard surface such as tile, wood, or linoleum.**
 - **This and most bath scales will not perform properly on thick-carpeted surfaces.**
 - **Do not get the scale wet – this can damage the scale and cause an unsafe condition.**

**Disclaimer: This scale tells jokes – it is not meant to offend anyone – especially you 😊
We hope you will enjoy this product for many years of fun weighing!**