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

VT-1967 W

Electronic body fat
analyzer scale

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PERSONAL SCALE VT-1967 W

The scale is intended to measure weight, calculate the body weight index and to measure fat, water and bone tissue percentage in a human body.

This scale uses bioelectric body impedance analysis «BIA» that allows to precisely calculate the body weight index and percentage of fat and bone tissues and water. The calculations are based on personal data (height, age, sex) input into the unit memory.

DESCRIPTION

1. Weighing platform
2. Contact areas
3. Display
4. Data selection and confirmation button «SET»
5. Data selection buttons «▲/▼»
6. Battery compartment lid
7. Weight measurement units selection button «kg/lb»

RECOMMENDATIONS FOR CUSTOMER ATTENTION!

Before using the unit, read these instructions carefully and keep them for future reference.

Use the unit according to its intended purpose only, as it is stated in this user manual. Mishandling the unit can lead to its breakage and cause harm to the user or damage to his/her property.

PRECAUTIONS

- **Using the excess weight analysis function is TOTALLY contraindicated to pregnant women.**
- **Using the excess weight analysis function is TOTALLY contraindicated to people having pacemakers and other medical devices with built-in batteries.**
- Handle your scale with care as it is a measuring instrument, do not subject it to high or low temperatures, high humidity, never expose it to direct sunlight and do not drop it.
- Use the scale away from heating units.
- Avoid getting of liquids into the scale body, this appliance is not waterproof. Do not store or use the scale in places with high humidity (above 80%), make sure that the scale does not come into contact with water or other liquids as this may cause failure of scale indication or its damage.

- Place the scale on a flat, dry and non-slip surface.
- The scale readouts will be more accurate if you do not use the scale on carpets or rugs.
- Do not step on the damp surface of the scale, do not step on the scale with wet feet, otherwise you can slip and get injured.
- Be careful while stepping on the scale: step on the scale first with one foot; make sure that the scale is steady and only then step on it with both feet.
- Step on the scale and dispose your weight evenly. Stand still while weighing. Do not jump on the scale surface.
- Replace the batteries in time. Strictly follow the polarity according to the markings when replacing the batteries.
- If you are not planning to use the scale for a long time, remove the batteries from the battery compartment.
- Do not leave children unattended to prevent using the unit as a toy.
- For children safety reasons, do not leave polyethylene bags used as packaging unattended.

Attention! Do not allow children to play with polyethylene bags or packaging film. Danger of suffocation!

- This unit is not intended for usage by children under 8 years of age.
- This unit is not intended for usage by people (including children over 8) with physical, neural and mental disorders or with insufficient experience or knowledge. Such persons may use this unit only if they are under supervision of a person who is responsible for their safety and if they are given all the necessary and understandable instructions concerning the safe usage of the unit and information about danger that can be caused by its improper usage.
- Do not repair the scale by yourself. Refer to authorized service centers for maintenance.
- To avoid damages, transport the unit in original package only.
- Keep the unit out of reach of children and disabled persons.

THE UNIT IS INTENDED FOR HOUSEHOLD USAGE ONLY

USING THE SCALE

After unit transportation or storage under cold (winter) conditions, it is necessary to keep it for

ENGLISH

at least two hours at room temperature before switching on.

- Unpack the unit, remove any stickers that can prevent unit operation.
- Clean the weighing platform (1) and the scale body with a soft slightly damp cloth, and then wipe dry.

Battery installation

- Remove the battery compartment lid (6), having preliminary unscrewed the fixing screw, and insert 2 «CR2032» batteries (included), following the polarity.
- Install the battery compartment lid (6) back to its place and tighten up the screw.

Note: If there is an isolation insert in the battery compartment, open the battery compartment lid (6), remove the isolation insert and install the battery compartment lid (6) back to its place.

Battery replacement

- The symbol «Lo» indicates that the battery is low, in this case replace the battery.
- Open the battery compartment lid (6), having preliminary unscrewed the fixing screw, replace the batteries with new CR2032 batteries, strictly following the polarity, close the battery compartment lid (6) and tighten up the screw.
- If you are not going to use the scale for a long time, remove the batteries from the battery compartment.

WHAT IS EXCESS WEIGHT?

It is known that excessive weight can negatively affect health. Excessive weight can cause obesity, heart problems, atherosclerosis and osteoporosis. Short weight may also lead to osteoporosis.

Osteoporosis is a systemic skeleton disease characterized by reduction of bone mass per volume unit and impairment of bone tissue structure, leading to increased bone fragility and risk of fracture by minimal trauma or even without it. Physical exercises and a moderate diet help to lose excess weight and this scale will help to control your weight.

It is known that in case of normal weight muscles and body tissues contain more water than fat, and that electric current passes through water better than through fat. Sending weak current through the body, the scale quickly and easily calculates the percentage of fat and water in a human body.

RECOMMENDATIONS

- Weigh yourself at one and the same time of the day.
- We recommend measuring the percentage of fat and water in the evening before evening meal.
- 2-3 hours should pass after the last meal.
- 10-12 hours should pass after intensive physical exercises.
- Percentage of fat in the body depends on different factors and can change during the day.
- When calculating percentage of fat and water of children younger than 6 years of age, the results may be inaccurate.
- Professional athletes, people inclined to edemas and people with nonstandard body proportions can show inaccurate results of calculation.

Note: *The percentage of fat and bone tissues and water and the body weight index can be calculated only if personal data (sex, age, height) are entered correctly and if you step on the scale barefoot.*

SWITCHING ON

- To switch the scale on, press on the platform (1), the symbols «0.0» and weight measurement units «kg or lb» will appear on the display (3).
- To select weight measurement units, press the «kg/lb» button (7), the following symbols will be shown on the display (3):
 - «kg» – weight in kilograms;
 - «lb» – weight in pounds;
- Place the scale on a flat solid surface; if the floor is uneven or covered with a carpet or rug, the weighing accuracy decreases.

WEIGHT MEASUREMENT

- Step on the weighing platform (1). Be careful while stepping on the scale: step on the scale first with one foot; make sure that the scale is steady and then step on it with both feet.
- The symbols «0.0 kg» or «0.0 lb» will appear on the display, depending on the selected measurement units.
- Keep your feet parallel to each other and distribute your weight evenly on the platform (1).
- Stay still while weighing. The weight indications on the display (3) will flash during

weighing, when the weight is fixed, the weight indications will glow constantly.

Attention!

- If the symbol «Err» appears on the display (3), it means that the scale is overloaded. Step off the scale immediately to avoid its damage.
- The scale will be switched off automatically in about 15 seconds after weight measurement.

ENTERING PERSONAL DATA

- Press the «SET» button (4), the number of one of the users «P-0 user», «P-1 user»..., «P-9 user», his/her height and the symbols «0.0 kg» or «0.0 lb» – depending on the selected measurement units – will appear on the display (3).
- Select the required user number by repeatedly pressing the buttons (5) «▲/▼». You can enter data for 10 users marked by symbols «P-0 user», «P-1 user»..., «P-9 user». Confirm the user selection by pressing the «SET» button (4), the flashing symbols of the user's number «P-0 user», «P-1 user»..., «P-9 user» and symbols «0.0 kg» or «0.0 lb» – depending on the selected measurement units – will appear on the display (3).
- Press the «SET» button (4) again, the flashing digital values of age and the symbols «age» will appear on the display (3). Set the user's age by repeatedly pressing or pressing and holding down one of the buttons (5) «▲/▼», the age range is from 10 to 100 years with the step 1 year. To save the user's age, press the «SET» button (4).
- Then the flashing symbol of the sex will be shown on display (3). Select the symbol  for women or  for men, pressing one of the buttons (5) «▲/▼». Confirm your selection by pressing the «SET» button (4).
- Then flashing digital values of height will be shown on the display (3), the digital values of height are entered in values corresponding to the selected weight measurement units. If you select the weight measurement units «kg», the user's height values are entered in «centimeters», possible height values are from 100 to 220 cm with 1 cm step. Set the height values by repeatedly pressing or pressing and holding down one of the buttons (5) «▲/▼», confirm your selection by pressing the «SET» button (4).

CALCULATION OF FAT, BONE AND MUSCLE TISSUE AND WATER PERCENTAGE IN A HUMAN BODY

- Press the «SET» button (4).
- By consequently pressing the buttons (5) «▲/▼», select the user number with your data «P-0 user», «P-1 user»..., «P-9 user».
- You can view the personal information by repeatedly pressing the «SET» button (4).
- After viewing your personal data, press the «SET» button (4) again, the symbols «0.0 kg» or «0.0 lb» will appear on the display (3), depending on the selected measurement units.

Note: If you do not need to view your personal data, then after having selected the user's number, wait till the symbols «0.0 kg» or «0.0 lb» – depending on the selected measurement units – will appear on the display (3).

- Step with your bare feet on the scale so that your feet are on the contact areas (2), stay still during the weight measurement and wait till the weight indications appear.
- Remain standing on the scale, the descending symbols «oooo» will appear on the display (3).
- After calculation, the following data will be shown on the display (3):
 - the user's weight in selected weight units (kg/lb);
 - fat percentage (BF),
 - water percentage (BW),
 - body weight index (BMI),
 - bone tissue percentage (BONE).
- The calculation results will be shown on the display for 30 seconds and then the unit will be switched off automatically.

Notes:

- IT IS NECESSARY to enter your personal data into the scale memory before measuring fat, water percentage, etc.
- The data on the percentage of fat and water may be obtained only if you stand on the scale barefoot.
- If the symbol «Err» appears on the display (3), it means that the scale is overloaded. Step off the scale immediately to avoid its damage. Maximal capacity of the scale is 180 kg.
- The symbol «Err2» (in this model only Err appears) indicates that fat percentage in the body is less than 5% or exceeds 85%.

ENGLISH

Change of fat percentage in a human body during the day

The indications of fat percentage in a human body are affected by such factors as lifestyle, physical activity, health condition and water percentage in the body.

You should take into account the following factors:

- the fat percentage indications may be less if measured after taking a bath or a shower,
- the fat percentage indications may be higher if measured after having a meal or a drink,
- percentage of liquid in a human body may decrease and thus change the fat percentage indications during and after illness and after physical exercises.

Age	Woman (% of fat)				Man (% of fat)			
	Low	Medium	High	Obese	Low	Medium	High	Obese
10-39	14-20	21-25	26-31	32-38	11-17	18-23	24-29	30-36
40-55	15-21	22-26	27-32	33-38	12-18	19-24	25-30	31-37
56-100	16-22	23-27	28-33	34-38	13-19	20-25	26-31	32-38

Note: The data in this table is given for recommendation purposes only and can vary, depending on individual parameters of the user.

Water percentage in a human body

Water percentage in a human body is the

percentage of liquid in the human body to total weight. Water has primary importance for human tissues and organs, since 50-65% of our body consist of water. With this scale you can control the percentage of water in your body and correct this important index.

Normal percentage of water in the body is 45-60% for women and 50-65% for men.

Age	Woman (% of water)			Man (% of water)		
	Low	Medium	High	Low	Medium	High
10-15	<57	57-67	>67	<58	58-72	>72
16-30	<47	47-57	>57	<53	53-67	>67
31-60	<42	42-52	>52	<47	47-61	>61
61-100	<37	37-47	>47	<42	42-56	>56

Note: The data in this table is given for recommendation purposes only and can vary, depending on individual parameters of the user.

Insufficient weight	Normal weight	Excessive weight
<20	20-25	>25

Body weight index (BMI)

The body weight index is the value that allows to assess the degree of conformity of a person's weight with his/her height and thus indirectly tell whether the body weight is insufficient, normal or excessive.

Apply the body weight index with caution, only for a rough estimate – for example, an attempt to assess by its means the body build of professional athletes can lead to incorrect results (high index value in this case is accounted for by the athlete's well developed musculature).

The body weight index is calculated according to the formula:

$I = h/m^2$, where m is the body weight in kilograms, h is the height in meters, the index is measured in kg/m².

Note: The data in this table is given for recommendation purposes only and can vary, depending on individual parameters of the user.

What is bone mass?

Bone tissues are the base of human skeleton; they are solid and resistant to external impact, protecting internal organs; together with the muscle tissue, bone tissues provide the ability of a human to move.

Bone mass measurement is based on calculating the minerals mass in a human body (calcium or other matters). Research has showed that development of muscle tissues makes the skeleton stronger. This fact should be taken into account when making up a diet and calculating physical activity to maintain the bones strong and healthy.

The bone mass practically doesn't change; however, some decrease is possible under many factors (weight, height, age, sex). That's why this manual does not contain standard values and recommendations concerning the bone tissue percentage.

Analysis of weighing results

As a rule, short-term changes of weight within several days are the result of minor loss of liquid in your body tissues. Analyze the results of regular weight measurements taken throughout a long time. Evaluation of the results is based on changes in the total weight, fat and water percentage in body tissues and depends on the duration of the period when the measurements are made.

Short-term indication changes during a day can be significantly different from the indication changes within a week or a month.

As a rule, short-term indication changes are connected with percentage of liquid in body tissues, and medium and long-term changes are caused by changes in the amount of fat and muscle bulk.

- For instance, short-term weight loss with simultaneously increased or constant fat percentage in body tissues may be caused by decreased liquid amount as a result of physical activity, after visiting sauna or by quick weight loss.
- Medium-term increase of weight, decrease of fat amount or constant fat amount in body tissues may be caused by increased muscle bulk.
- If you notice simultaneous decrease of total weight and fat percentage, it means that your diet works, i.e. you lose the fat content in your body tissues.
- The best results are achieved if you combine diet with various physical exercises. In this case you can observe the medium-term increase of muscle mass in body tissues.

CLEANING AND CARE

1. Wipe the scale with a soft, slightly damp cloth and then wipe it dry.
2. Never use abrasives or solvents.
3. Never immerse the scale into water or other liquids.
4. Remove the batteries from the battery compartment before taking the unit away for storage.
5. Keep the scale away from children in a dry cool place.

DELIVERY SET

1. Scale – 1 pc.
2. «CR2032» battery – 2 pc.
3. Instruction manual – 1 pc.

TECHNICAL SPECIFICATIONS

Power supply: 2 CR2032 batteries, 3 V

Max. weight: 200 kg.

The manufacturer preserves the right to change the specifications of the unit without a preliminary notification.

Unit operating life is 3 years

Guarantee

Details regarding guarantee conditions can be obtained from the dealer from whom the appliance was purchased. The bill of sale or receipt must be produced when making any claim under the terms of this guarantee.

This product conforms to the EMC-Requirements as laid down by the Council Directive 2004/108/EC and to the Low Voltage Regulation (2006/95/EC)

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© ООО ГОЛДЕР-ЭЛЕКТРОНИКС, 2014
© GOLDER-ELECTRONICS LLC, 2014